


Rectal Self-Swab Collection Instructions




A landmark in prevention



Step 1.

Open kit and remove tube and package with green writing. Remove the swab with the **BLUE** shaft. **USE BLUE SHAFT SWAB ONLY.**



Step 3.



Remove cap from test tube. Place swab in test tube. Do not puncture the foil cap.

Break swab shaft at the score mark.

Step 2.

Insert swab 1 inch into the anus and turn for 5 – 10 seconds.

If needed, before inserting swab, wet swab with water or saline solution.



Step 4.

Put cap back tightly on test tube to prevent any leaking. Try not to splash the liquid out the tube.

Step 5.

Discard wrapper and unused swab. **Wash your hands.** Return the tube to the health worker.